

## MolokClassic BioSystem 500 C

Used in 800 l MolokClassic® container to replace the lifting bag when collecting organic waste. Can be used for other waste types, such as metal and glass as well.

Well body material HD Polyethylene

Volume 447 l

+ liquid container 48 l

SWL 400 kg Total permissible mass 430 kg

Rotation molded lifting container made from HD polyethylene. Plastic locking mechanism and hinges from durable textile material. Midsole between bottom cup (48 l) and container body (447 l) separates liquids from the solid waste.

## Additional information

- Organic waste bag can be used inside the container to keep it clean. The disposable organic waste bag is fastened with plastic fixing cone.
- Emptying is done by lifting loops. The container is emptied by opening the locking device, which opens up the bottom cup and makes the waste drop to the the emptying truck.

